GRIND GUARD

Anti Teeth Grinding Mouthpiece

Scan for
Instructions

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EN

PRODUCT DESCRIPTION

The Grind Guard is a personal-use night guard designed to help protect teeth from the effects of nighttime grinding (bruxism). When softened in warm water, the thermoplastic material molds to the user's dental shape, providing a cushioning layer between the upper and lower teeth. This product is not a medical device and is not intended to diagnose or treat any dental condition. It is suitable for light-to-moderate grinding support and routine nighttime use when properly fitted and cleaned.

INTENDED USE

- Helps reduce wear on teeth caused by nighttime grinding.
- Provides a protective barrier between upper and lower teeth.
- Helps relieve mild tension from jaw clenching during sleep.

This product is not intended for individuals with severe bruxism or TMJ disorders unless advised by a dentist or healthcare professional.

INSTRUCTIONS FOR USE

Read all instructions and warnings carefully.

Fitting the Grind Guard

- Place the grind guard in a heat-resistant container suitable for immersion in water.
- Fill the container with hot water (70–80 °C) and fully submerge the grind guard. Let it sit for 15–25 seconds until the material begins to soften.
- Carefully lift the grind guard from the water using wooden tongs, chopsticks, or another heat-safe utensil to avoid contact burns.
- Shake off excess water and place the grind guard on a paper towel for 3–5 seconds to cool slightly. Insert the grind guard into your mouth over your upper teeth.
- Bite down firmly, molding the guard to the contours of your teeth and gums. Press with your fingers along the front and back teeth to secure a tight fit.
- After creating a clear dental impression, immerse the grind guard in cold water for 30 seconds to set the mold.

If the fit is not correct: Reheat and repeat the process. The device can be refitted several times without affecting performance.

Using the Device

- Insert the Grind Guard before sleep.
- Ensure it sits securely and does not cause pressure or irritation.
- Remove immediately if discomfort or pain occurs.

Cleaning & Care

Rinse with warm water after each use. Do not use boiling water. Store the guard in the included ventilated case. Once a week, clean using mild soap water. Do not use bleach, alcohol, or abrasive cleaners.

▲ SAFETY INFORMATION & WARNINGS

- Not recommended for children under 12 years of age.
- Do not use if you wear braces, have loose or shifting teeth, or have dental implants that are still healing.
- Not suitable for individuals with severe bruxism, jaw disorders, or TMJ pain. Consult a dentist.
- Discontinue use if you experience:
 - Tooth or jaw pain
 - Gum irritation
 - Excessive drooling
 - Worsening sleep quality
- Keep out of reach of children when not in use.
- Do not use if the device shows visible damage, tears, or deformation
- Do not share the device with another person.
- Consult a dentist if discomfort or pain persists.

MATERIALS

- Food Grade EVA (BPA-free, latex-free)
- Storage case: Polypropylene (PP)

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